



## NEED A LIFT?

Golden lift chairs offer the utmost in quality, comfort and safety for you or your loved one. Our lift chairs are constructed by hand in Old Forge, Pennsylvania with a solid hardwood maple frame with the best warranty in the industry. Our chairs are packed with patented technology and Golden offers the widest selection of models, sizes, fabrics and features on the market. Visit a Golden dealer to be sized and customize your comfort today!



CONTACT US



[www.goldentechnologies.ca](http://www.goldentechnologies.ca)



# LIFT CHAIR SIZING GUIDE





## THE PERFECT FIT


Find your perfect fit in a Golden power lift recline chair. Start the sizing process by removing your coat, shoes and belongings and relax as if you are at home. Sit in a chair with straight posture and your feet firmly planted on the floor. Complete measurements A-F and have a Golden lift chair expert help you find the perfect fit for Customized Comfort.

### Giving a Lift Chair as a Gift?

Measure the main users favourite chair and we can help you find the closest match!



## MEASURING FOR CUSTOM COMFORT



### A - SEAT TO TOP OF BACK

Measure from the base of the seat to the top of the head in inches

### B - SEAT DEPTH

Measure from the back of the seat to the back of the knee in inches

### C - FLOOR TO TOP OF SEAT

Measure from the floor to the top of the back of the knee in inches

### D - WIDTH BETWEEN ARMS

Measure from the left hip to the right hip while seated and add 2 inches

### E - HEIGHT

The height of the main user of the lift chair in feet

### F - WEIGHT

The weight of the main user of the lift chair in pounds

## IMPORTANCE OF SIZING

### COMFORT

Finding the perfect fit for your new lift chair will provide the most comfortable option available. Golden offers the most sizes in the industry to customize your comfort!

### SAFETY

Your safety is our top priority! If a lift chair is too big or too small it can become a safety hazard when transferring in and out of the chair. We want to make sure you find the perfect fit for safe transfers everytime you enter and exit your chair.

### LONGEVITY

Did you know that if a lift chair is too big or too small it can cause mechanical issues over time? We want to help you find the perfect fit to encourage comfort and not headaches due to mechanical breakdown.